



# Sepor Vibra Pad

## CONTACT US

Phone: 310-830-6601

FAX: 310-830-9336

Email: [info@sepor.com](mailto:info@sepor.com) &  
[services@sepor.com](mailto:services@sepor.com)

## LOCATIONS

### WEST COAST

718 N. Fries Avenue  
Wilmington, CA 90744

### EAST COAST

3740 NW 124th Ave  
Coral Springs, FL 33065



## Description:

A specially designed plastic-coated pad mounted on a variable speed vibrator to perform many laboratory operations. Originally designed to accommodate the SEPOR Micro splitter and from one to six 3" diameter nested testing sieves, this unit has found many applications such as compaction, settling, bubble removal, etc.

For 115v/50-60hz. Variable control, 6-ft cord and plug, resilient mounting.

The Vibra Pad is 5 1/4" square with 4 pillars, two anchored studs, wing nuts, washers, and garter springs.