



Gyratory Sieve Shaker

CONTACT US

Phone: 310-830-6601

FAX: 310-830-9336

Email: info@sepor.com
& services@sepor.com


LOCATIONS

WEST COAST

718 N. Fries Avenue
Wilmington, CA 90744

EAST COAST

3740 NW 124th Ave
Coral Springs, FL 33065



Wet Sieve Kit Drain Pan



Wet Kit Spray Cover

Wet Sieve Kit Spray Cover

Description:

The Gyratory Sifters are designed for laboratory or quality control testing, where accurate and reproducible results are important. The 8" Gyratory Sifter uses standard 8" lab sieves and has a capacity for 8 full height or 16 half height sieves. These lab sifters must be suspended from either a ceiling or a stable overhead support. They are popular for wet sieving, which requires the purchase of the optional wet sieve kit, pictured above.

Construction:

The 8" Gyratory Sifter has a 1/6 HP/110 V/1 Ph/60 Hz or 1/6 HP/220 V/1 Ph/50 Hz motor and drive.

Catalog Number	Description	Ship Weight
040B-021	8" Gyratory Sieve Shaker	55 Lbs
040B-023	Wet Sieve Kit For 8" Sieve Shaker	15 Lbs.

This gyratory sieve shaker is primarily used for wet sieving, and most labs will place a 5-gallon bucket beneath the drain pan, with rubber tubing attached, to capture the fines and water. The fines can be de-watered in **a Sepor pressure filter, which come in different sizes depending upon your needs.**