



## Lab Rare Earth Roll Magnetic Separator L/P 10-30

### CONTACT US

Phone: 310-830-6601

FAX: 310-830-9336

Email: [info@sepor.com](mailto:info@sepor.com)

& [services@sepor.com](mailto:services@sepor.com)


### LOCATIONS

#### WEST COAST

718 N. Fries Avenue  
Wilmington, CA 90744

#### EAST COAST

3740 NW 124th Ave  
Coral Springs, FL 33065



### Description:

The Laboratory High-Force Magnetic Separator L/P 10-30 is designed to separate paramagnetic materials from non-magnetic materials. The magnetic zone is 200mm wide and has a dual magnetic roll configuration of 4:1.3/ 4:4. The magnetic roll is made with the highest strength Neodymium magnets. It includes the hopper system, two splitters and the frame, all made of stainless steel.

This separator can operate continuously and has a variable speed drive system (please specify electrical requirements). It includes closed loop speed control to maintain roll speed at desired setting and a digital RPM read-out. The vibratory feeder is mounted in the frame along with the unit's feeder. Three stainless steel collection pans and three belts are included as standard.

Catalog Number	Description	Power	Ship Weight	Dimensions (in.)
<b>060C-010</b>	Laboratory Rare Earth Roll Magnetic Separator L/P 10-30	220v/1/50-60Hz	475	41 x 38 x 42